

NKX1-2 antibody



Catalog Number: 113202

Product name

NKX1-2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NKX1-2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse brain tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001146340). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113202(NKX1-2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:113202(NKX1-2 Antibody) at dilution of 1:50 (under 40x lens)