

# PRKG1 antibody

Catalog Number: 114202

## Product name

PRKG1 antibody

## Specificity

Human, Rat, Mouse; other species not tested.

## Antibody description

PRKG1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in rat brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 65-78 kDa

## Preparation

This antibody was obtained by immunization of PRKG1 recombinant protein (Accession Number: NM\_006258). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114202(PRKG1 Antibody) at dilution of 1:600



IP Result of anti-PRKG1 (IP:Catalog No:114202, 4ug; Detection:Catalog No:114202 1:600) with mouse brain tissue lysate 4000ug.