

SERPINB13 antibody



Catalog Number: 115202

Product name

SERPINB13 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SERPINB13 Rabbit Polyclonal antibody. Positive IHC detected in human skin. Positive IF detected in HepG2.

Preparation

This antibody was obtained by immunization of SERPINB13 recombinant protein (Accession Number: NM_012397). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human skin using Catalog No:115202(SERPINB13 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human skin using Catalog No:115202(SERPINB13 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using SERPINB13 antibody Catalog No:115202 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).