

# FAK antibody

Catalog Number: 107202

## Product name

FAK antibody

## Specificity

Human; other species not tested.

## Antibody description

FAK Mouse Monoclonal antibody. Positive IF detected in A549 cells. Positive WB detected in MCF-7 cells, human testis tissue. Observed molecular weight by Western-blot: 110 kDa

## Preparation

This antibody was obtained by immunization of FAK recombinant protein (Accession Number: BC028733). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2c

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with (FAK Antibody) at dilution of



1:1000

Immunofluorescent analysis of (10% Formaldehyde) fixed A549 cells using Catalog No:107202(FAK Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)