

Anti-beta Catenin (7E2) Mouse antibody



Catalog Number: 168202

Product name

WB: 1/1000-2000

Anti-beta Catenin (7E2) Mouse antibody

IHC: 1/200-500

Specificity

IF: 1/200

Human, Mouse, Rat

Antibody description

beta Catenin (7E2) Mouse monoclonal antibody

Preparation

Antigen: Recombinant Protein of Catenin- β

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


WB, ICC/IF, IHC-P

Dilutions

Validations



Western blot analysis of 1) HeLa, 2) 293T, 3) Mouse Liver Tissue, 4) Rat Liver Tissue, with β -Catenin Mouse mAb diluted at 1:2,000.



Western blot detection of β -Catenin in T47D, 3T3, COS7 and C6 cell lysates using β -Catenin mouse mAb (1:3000 diluted). Predicted band size: 92KDa. Observed band size: 92KDa.