

FAM84B antibody

Catalog Number: 110602

Product name

FAM84B antibody

Specificity

Human; other species not tested.

Antibody description

FAM84B Rabbit Polyclonal antibody. Positive IP detected in HUVEC cells. Positive WB detected in MDA-MB-453s cells, Jurkat cells. Positive IHC detected in human liver cancer tissue, human heart tissue, human kidney tissue, human liver tissue, human lung tissue, human oesophagus cancer tissue, human ovary tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of FAM84B recombinant protein (Accession Number: XM_017013107). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations




MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:110602(FAM84B Antibody) at dilution of 1:600




IP Result of anti-FAM84B (IP:Catalog No:110602, 4ug; Detection:Catalog No:110602 1:300) with HUVEC cells lysate 3600ug.

FAM84B antibody

Catalog Number: 110602



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:110602(FAM84B Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:110602(FAM84B Antibody) at dilution of 1:200 (under 40x lens).