

Mouse IGFBP-2 / IGFBP2 (Fc Tag) recombinant protein



Catalog Number: 501602

General Information

Protein Construction

A DNA sequence encoding the mouse Igfbp2 (NP_032368.2) (Met1-Gln305) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> (65.9+24.5) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 35

Molecular Mass

The recombinant mouse Igfbp2 consists 512

amino acids and predicts a molecular mass of 56.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IGFBP-2 / IGFBP2 Protein (Fc Tag)