

PATZ1 antibody



Catalog Number: 113602

Product name

PATZ1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PATZ1 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, Raji cells. Observed molecular weight by Western-blot: 70-74 kDa

Preparation

This antibody was obtained by immunization of PATZ1 recombinant protein (Accession Number: NM_032051). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113602(PATZ1 Antibody) at dilution of 1:1000