

SRPX2 antibody



Catalog Number: 107602

Product name

SRPX2 antibody

Specificity

Human, Pig; other species not tested.

Antibody description

SRPX2 Mouse Monoclonal antibody. Positive WB detected in pig lung tissue, A549 cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of SRPX2 recombinant protein (Accession Number: NM_014467). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



pig lung tissue were subjected to SDS PAGE followed by western blot with (SRPX2 Antibody) at dilution of 1:1000