

# Anti-GNA13 Rabbit antibody



Catalog Number: 168602

## Product name

Anti-GNA13 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GNA13 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human G protein alpha 13

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000-1/5000

IHC: 1/200

## Validations



Western blot detection of GNA13 [KO Validated] in HeLa, Jurkat cell lysates using GNA13 Rabbit mAb [KO Validated](1:500 diluted). Predicted band size: 44kDa. Observed band size: 44kDa.