

Human UPRT (His tag) recombinant protein



Catalog Number: 518602

General Information

Gene Name Synonym

Uracil phosphoribosyltransferase homolog, UPP, FUR1, UPRT

Protein Construction

A DNA sequence encoding the Human UPRT Met1-Asp309 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation

Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl pH 8.0.


Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE

kDa MK R NR



Human UPRT (His tag) recombinant protein