

HRP-conjugated Mouse anti GST-Tag mAb



Catalog Number: 166302

Product name

HRP-conjugated Mouse anti GST-Tag mAb

Antibody description

Monoclonal antibody to HRP-conjugated Mouse anti GST-Tag

Preparation

Antigen: Recombinant protein of GST-Tag

Formulation

PBS, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB 1:2000 - 1:10000

Validations



Western blot - HRP-conjugated Mouse anti GST-Tag mAb

Western blot analysis of recombinant GST protein with HRP-conjugated Mouse anti GST-Tag antibody at dilution of 1:5000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.