

Human MAN1B1/UNQ747/PRO1477 (His tag) recombinant protein



Catalog Number: 519302

General Information

Gene Name Synonym

Endoplasmic Reticulum Mannosyl-Oligosaccharide 1,2-Alpha-Mannosidase, ER Alpha-1,2-Mannosidase, ER Mannosidase 1, ERMan1, Man9GlcNAc2-Specific-Processing Alpha-Mannosidase, Mannosidase Alpha Class 1B Member 1, MAN1B1

Protein Construction

A DNA sequence encoding the Human MAN1B1/UNQ747/PRO1477 Asp106-Ala699 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 50mM TrisHCL, 10mM reduced Glutathione, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human MAN1B1/UNQ747/PRO1477 (His tag)
recombinant protein