

# RAD54L2 Polyclonal Antibody



Catalog Number: 161702

---

## Product name

RAD54L2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to RAD54L2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1178-1467 of human RAD54L2 (NP\_055921.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - RAD54L2 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RAD54L2 antibody . Blue: DAPI for nuclear staining.