

Anti-SPINK4 antibody



Catalog Number: 102702

Product name

Anti-SPINK4 antibody

Immunogen

[Human SPINK4 \(His Tag\) recombinant protein](#)

Specificity

Human SPINK4

Antibody description

Rabbit polyclonal to SPINK4

Preparation

Produced in rabbits immunized with purified, recombinant Human SPINK4 (rh SPINK4; NP_055286.1; Met 1-Cys 86). SPINK4 specific IgG was purified by human SPINK4 affinity chromatography

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human SPINK4. The detection limit for Human SPINK4 is approximately 0.00975 ng/well.


IHC-P: 0.1-2 µg/mL

Validations



SPINK4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human SPINK4 in human pancreas with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded



sections).

SPINK4 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human SPINK4 in

Anti-SPINK4 antibody



Catalog Number: 102702

human skin with rabbit polyclonal antibody (1 μ g/mL, formalin-fixed paraffin embedded sections).