

Wnt5a Polyclonal Antibody



Catalog Number: 163702

Product name

Wnt5a Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to Wnt5a

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 101-380 of human Wnt5a (NP_003383.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - Wnt5a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using Wnt5a antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - Wnt5a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using Wnt5a antibody at dilution of 1:100 (40x lens).