

PCDHGC3 antibody



Catalog Number: 113702

Product name

PCDHGC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCDHGC3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue.

Preparation

This antibody was obtained by immunization of PCDHGC3 recombinant protein (Accession Number: NM_002588). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113702(PCDHGC3 antibody) at dilution of 1:100 (under 10x lens)