

# Anti-p23/PTGES3 Rabbit antibody



Catalog Number: 167702

## Product name

Anti-p23/PTGES3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

p23/PTGES3 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human p23/PTGES3

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF, IP

## Dilutions

WB: 1/2000-1/100001/20

IP: 1/20

## Validations



Western blot detection of p23/PTGES3 in HeLa,A549,HL-60 using p23/PTGES3 Rabbit mAb(1:1000 diluted)