

# C20orf4 antibody



Catalog Number: 108702

## Product name

C20orf4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

C20orf4 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 40-45 kDa

## Preparation

This antibody was obtained by immunization of C20orf4 recombinant protein (Accession Number: NM\_015511). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108702(C20orf4 Antibody) at dilution of 1:300