

# DR4 antibody

Catalog Number: 110012

## Product name

DR4 antibody

## Specificity

Human; other species not tested.

## Antibody description

DR4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 44-50 kDa

## Preparation

This antibody was obtained by immunization of DR4 recombinant protein (Accession Number: NM\_003844). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110012(DR4 antibody) at dilution of 1:1000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110012(DR4 Antibody) at dilution of 1:20000