

NAB2 antibody

Catalog Number: 113012

Product name

NAB2 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

NAB2 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human brain tissue, K-562 cells, mouse thymus tissue. Positive IP detected in K-562 cells. Positive IHC detected in human brain tissue, human kidney tissue. Observed molecular weight by Western-blot: 57 kDa

Validations




Preparation

This antibody was obtained by immunization of NAB2 recombinant protein (Accession Number: NM_005967). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113012(NAB2 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:113012(NAB2 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded

NAB2 antibody

Catalog Number: 113012

human brain using Catalog No:113012(NAB2 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-NAB2 (IP:Catalog No:113012, 4ug; Detection:Catalog No:113012 1:700) with K-562 cells lysate 3200ug.

