

# Anti-Rxrp3 antibody

Catalog Number: 175012

---

## Product name

Anti-Rxrp3 antibody

## Specificity

Mouse

## Antibody description

Rabbit polyclonal antibody to Rxrp3

## Preparation

This antigen of this antibody was synthetic peptide within mouse relaxin 3 receptor 1 aa 300-352.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC, IHC-P, FC

## Dilutions

WB: 1:5,000~1:10,000

ICC: 1:100~1:200

## Validations




Fig1: Western blot analysis on mouse spinal cord (A) and mouse brain (B) using anti- Relaxin 3 receptor 1 polyclonal antibody.