

Anti-CD32a / Fc gamma RIIA / FCGR2A antibody



Catalog Number: 106012

Product name

Anti-CD32a / Fc gamma RIIA / FCGR2A antibody

Immunogen

[Human CD32a / Fc gamma RIIA / FCGR2A \(167 Arg, His Tag\) recombinant protein](#)

Specificity

Human CD32a / FCGR2A

Antibody description

Rabbit monoclonal to CD32a / Fc gamma RIIA / FCGR2A

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD32a / FCGR2A (rh CD32a / FCGR2A; AAA35827.1; Met1-Ile218).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 20-50 µg/mL

Validations



Items	Lanes	A
Sample (whole cell lysate)		He-1a
Sample Volume (µg/lane)		60
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		20-50 µg/ml
Secondary Antibody		DyLight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size: 36 kDa Observed band size: 37 kDa Additional bands at: 53 kDa (We are unsure as to the identity of this extra band. It may be a non-specific banding.)

CD32a / FCGR2A Antibody, Rabbit MAb, Western blot