

Anti-MAPK6/ERK-3 Rabbit antibody



Catalog Number: 167012

Product name

Anti-MAPK6/ERK-3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MAPK6/ERK-3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human MAPK6

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000-1/5000

IHC: 1/20

FC: 1/20-1/100

IP: 1/20

Validations



Western blot detection of MAPK6 in 3T3,HeLa cell lysates using MAPK6 Rabbit pAb(1:1000 diluted).Predicted band size:83kDa.Observed band size:105kDa.