

# Retifanlimab reference antibody



Catalog Number: 177012

## Product name

Retifanlimab reference antibody

## Specificity

Human PDCD1 / PD-1 / CD279

## Antibody description

Anti-PDCD1 / PD-1 / CD279 Antibody  
(Retifanlimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1×PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-PDCD1 / PD-1 / CD279 Antibody (retifanlimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-PDCD1 / PD-1 / CD279 Antibody (retifanlimab) is more than 95%, determined by SEC-HPLC.