

ZNF611 antibody



Catalog Number: 117012

Product name

ZNF611 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF611 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 81kd

Preparation

This antibody was obtained by immunization of ZNF611 recombinant protein (Accession Number: BC000918). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117012(ZNF611 antibody) at dilution of 1:500