

CD24 antibody



Catalog Number: 109012

Product name

CD24 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD24 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, Jurkat cells, mouse spleen tissue, Raji cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of CD24 recombinant protein (Accession Number: BC007674). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109012(CD24 antibody) at dilution of 1:500