

# Anti-SerpinB4 / SCCA2 antibody



Catalog Number: 100412

## Product name

Anti-SerpinB4 / SCCA2 antibody

## Immunogen

[Human SerpinB4 / SCCA2 \(His Tag\) recombinant protein](#)

## Specificity

Human SerpinB4 / SCCA2

## Antibody description

Rabbit monoclonal to SerpinB4 / SCCA2

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human SerpinB4 / SCCA2 (rh SerpinB4 / SCCA2; P48594; Met1-Pro390).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 1-10 µg/mL

## Validations



SerpinB4 / SCCA2 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human SERPINB4 in human skin with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

SerpinB4 / SCCA2 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human SERPINB4 in human esophageal carcinoma with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).