

CTRB1 Polyclonal Antibody



Catalog Number: 162412

Product name

CTRB1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to CTRB1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 79-263 of human CTRB1 (NP_001897.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - CTRB1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CTRB1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunohistochemistry - CTRB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using CTRB1 antibody at dilution of



1:100 (40x lens).

CTRB1 Polyclonal Antibody



Catalog Number: 162412

Immunohistochemistry - CTRB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using CTRB1 antibody at dilution of 1:100 (40x lens).