

ARPC1A Polyclonal Antibody



Catalog Number: 163412

Product name

ARPC1A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ARPC1A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 171-350 of human ARPC1A (NP_006400.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Immunohistochemistry - ARPC1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using ARPC1A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ARPC1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using ARPC1A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ARPC1A Polyclonal


ARPC1A Polyclonal Antibody



Catalog Number: 163412

Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using ARPC1A antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - ARPC1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using ARPC1A antibody at dilution of 1:100 (40x lens).