

# TRPM6 antibody

Catalog Number: 116412

---

## Product name

TRPM6 antibody

## Specificity

Human; other species not tested.

## Antibody description

TRPM6 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_017662).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of HEK-293 cells using Catalog No:116412(TRPM6 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG