

Anti-CTNNBL1 Mouse antibody



Catalog Number: 168412

Product name

Anti-CTNNBL1 Mouse antibody

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

Specificity

Human

ELISA: 1/10000

Antibody description

CTNNBL1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CTNNBL1 (AA: 390-557) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG2a

Applications

WB, FC


Dilutions

Validations



Flow cytometric analysis of HeLa cells using CTNNBL1 mouse mAb (green) and negative

control (red).



Western blot detection of CTNNBL1 in A549, MCF7 and HeLa cell lysates using CTNNBL1 mouse mAb (1:500 diluted). Predicted band size: 65KDa. Observed band size: 65KDa.