

GABRG1 antibody

Catalog Number: 110812

Product name

GABRG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GABRG1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of GABRG1 recombinant protein (Accession Number: NM_173536). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:110812(GABRG1 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110812(GABRG1 antibody) at dilution of 1:100 (under 40x lens)