

Human LAG3 recombinant protein (Native)



Catalog Number: 514812

General Information

Gene Name Synonym

LAG-3; Protein FDC; CD antigen CD223

Protein Construction

A DNA sequence encoding the human LAG3 (P18627) (Met1-Arg440) was expressed with an initial Met.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Val 29

Molecular Mass

The recombinant human LAG3 consists of 412 amino acids and predicts a molecular mass of 44.7 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



Human LAG3/LAG-3/CD223 Protein 15562