

Anti-IGFN1 antibody

Catalog Number: 175812

Product name

Anti-IGFN1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to IGFN1

Preparation

This antigen of this antibody was recombinant protein of igfn1

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:2000-1:5000

IHC-P: 1:200-1:500

Validations




Fig1:: Western blot analysis of 1) Pig Skeletal MuscleTissue, 2) Rat Skeletal Muscle Tissue, 3) Rat Skeletal Muscle Tissue using IGFN1 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded Mouse Skeletal MuscleTissue using IGFN1 Polyclonal Antibody.