

NELL1 antibody



Catalog Number: 113112

Product name

NELL1 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IHC: 1:20-1:200

[Human NELL1 Recombinant protein \(His tag\)](#)


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NELL1 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human kidney tissue, human liver tissue. Positive WB detected in mouse colon tissue, human testis tissue. Observed molecular weight by Western-blot: 110-120 kDa




Preparation

This antibody was obtained by immunization of NELL1 recombinant protein (Accession Number: NM_006157). Purification method: Antigen affinity purified.

mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:113112(NELL1 antibody) at dilution of 1:200

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113112(NELL1 Antibody) at dilution of 1:50 (under 10x lens)

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions



NELL1 antibody



Catalog Number: 113112

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog

No:113112(NELL1 Antibody) at dilution of 1:50 (under 40x lens)