

# NELL1 antibody



Catalog Number: 113112

## Product name

NELL1 antibody

## Immunogen

[Human NELL1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NELL1 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human kidney tissue, human liver tissue. Positive WB detected in mouse colon tissue, human testis tissue. Observed molecular weight by Western-blot: 110-120 kDa

## Preparation

This antibody was obtained by immunization of NELL1 recombinant protein (Accession Number: NM\_006157). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:113112(NELL1 antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113112(NELL1 Antibody) at dilution of 1:50 (under 10x lens)



# NELL1 antibody



Catalog Number: 113112

---

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog

No:113112(NELL1 Antibody) at dilution of 1:50 (under 40x lens)