

PPP3CC antibody



Catalog Number: 114112

Product name

PPP3CC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP3CC Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue, MCF7 cells, mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of PPP3CC recombinant protein (Accession Number: NM_005605). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114112(PPP3CC antibody) at dilution of 1:500