

TMEM194A antibody



Catalog Number: 116112

Product name

TMEM194A antibody

Clonality

Polyclonal

Immunogen

[Human TMEM194A Recombinant protein \(His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

TMEM194A Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, mouse liver tissue. Observed molecular weight by Western-blot: 40-42, 50-55 kDa


Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of TMEM194A recombinant protein (Accession Number: NM_015257). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116112(TMEM194A Antibody) at dilution of 1:1000