

ALK,CD246 antibody



Catalog Number: 107112

Product name

ALK,CD246 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

[Human ALK,CD246 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ALK,CD246 Mouse Monoclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of ALK,CD246 recombinant protein (Accession Number: NM_004304). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


Validations



Immunohistochemistry of paraffin-embedded human lymphoma slide using (ALK,CD246 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human lymphoma slide using (ALK,CD246



Antibody) at dilution of 1:50

Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog

ALK,CD246 antibody



Catalog Number: 107112

No:107572(ALK,CD246 Antibody) at dilution of 1:1000