

# ALK,CD246 antibody



Catalog Number: 107112

## Product name

ALK,CD246 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


## Immunogen

[Human ALK,CD246 Recombinant protein \(His tag\)](#)

## Validations

## Specificity

Human; other species not tested.




## Antibody description

ALK,CD246 Mouse Monoclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in Recombinant protein.

Immunohistochemistry of paraffin-embedded human lymphoma slide using (ALK,CD246 Antibody) at dilution of 1:50

## Preparation

This antibody was obtained by immunization of ALK,CD246 recombinant protein (Accession Number: NM\_004304). Purification method: Protein A purified.



Immunohistochemistry of paraffin-embedded human lymphoma slide using (ALK,CD246

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type


Mouse IgG2a

## Applications

ELISA, IHC, WB

## Dilutions

Recommended Dilution:



Antibody) at dilution of 1:50

Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog

# ALK,CD246 antibody



Catalog Number: 107112

---

No:107572(ALK,CD246 Antibody) at dilution of 1:1000