

Anti-B3GALT4 Mouse antibody



Catalog Number: 169112

Product name

Anti-B3GALT4 Mouse antibody

Specificity

Human

Antibody description

B3GALT4 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human B3GALT4 (AA: 191-359) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2a


Applications

FC


Validations



Western blot analysis using B3GALT4 mAb against HEK293 (1) and B3GALT4 (AA



Flow cytometric analysis of PANC-1 cells using B3GALT4 mouse mAb (green) and negative



control (red).

Western blot analysis using B3GALT4 mAb against

Anti-B3GALT4 Mouse antibody



Catalog Number: 169112

Western blot analysis using B3GALT4 mouse mAb against PANC-1 (1), PC-3 (2) cell lysate.

