

OTOP3 antibody



Catalog Number: 113512

Product name

OTOP3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

OTOP3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse brain tissue. Observed molecular weight by Western-blot: 90-100kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_178233). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113512(OTOP3 antibody) at dilution of 1:500