

WNT7A Polyclonal Antibody



Catalog Number: 164512

Product name

WNT7A Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to WNT7A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 32-170 of human WNT7A (NP_004616.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations



Western blot - WNT7A Polyclonal Antibody

Western blot analysis of extracts of mouse kidney, using WNT7A antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - WNT7A Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using WNT7A antibody.