

UBD antibody

Catalog Number: 116512

Product name

UBD antibody

Specificity

Human; other species not tested.

Antibody description

UBD Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue, Jurkat cells, Raji cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 18 kDa

Preparation

This antibody was obtained by immunization of UBD recombinant protein (Accession Number: NM_006398). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:116512(UBD antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116512(UBD antibody) at dilution of 1:100 (under 10x lens)