

Anti-Calbindin Rabbit antibody



Catalog Number: 167512

Product name

Anti-Calbindin Rabbit antibody

Specificity

Mouse, Rat

Antibody description

Calbindin Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse Calbindin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/200

IP: 1/20

Validations



Western blot detection of Calbindin in Rat Brain lysates using Calbindin Rabbit pAb(1:1000 diluted).Predicted band size:30kDa.Observed band size:30kDa.