

Anti-AMPK beta 1 Mouse antibody



Catalog Number: 168512

Product name

Anti-AMPK beta 1 Mouse antibody

Ig Type

Mouse IgG

Specificity

Human

Applications

WB

Antibody description

AMPK beta 1 Mouse monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: Recombinant Protein of AMPK β 1.

Validations

Formulation


Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 50% glycerol 0.1%BSA

Storage

Store at -20°C

Clonality

Monoclonal



Western blot detection of AMPK β 1 in HeLa, COS7 cell lysates using AMPK beta 1 Mouse mAb(1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.