

# Human PPIH/CYP20/CYPH (His tag) recombinant protein



Catalog Number: 519512

## General Information

### Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase H, PPIase H, Rotamase H, Small Nuclear Ribonucleoprotein Particle-Specific Cyclophilin H, CypH, U-snRNP-Associated Cyclophilin SnuCyp-20, USA-CYP, PPIH, CYP20, CYPH

### Protein Construction

A DNA sequence encoding the Human PPIH/CYP20/CYPH Met1-Met177 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human PPIH/CYP20/CYPH (His tag) recombinant protein