

# TCN1 antibody

Catalog Number: 115912

## Product name

TCN1 antibody

## Specificity

Human; other species not tested.

## Antibody description

TCN1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human tonsillitis tissue. Positive WB detected in human saliva tissue. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of TCN1 recombinant protein (Accession Number: BC018632). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



human saliva tissue were subjected to SDS PAGE followed by western blot with Catalog No:115912(TCN1 Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115912(TCN1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115912(TCN1 Antibody) at dilution of 1:50

# TCN1 antibody



Catalog Number: 115912

---

(under 40x lens)