

# Anti-Glutathione Reductase Rabbit antibody



Catalog Number: 166912

## Product name

Anti-Glutathione Reductase Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, FC

## Antibody description

Glutathione Reductase Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

## Preparation


Antigen: A synthetic peptide of human Glutathione Reductase

FC: 1/20

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of Glutathione Reductase in 3T3,HeLa cell lysates using Glutathione Reductase Rabbit pAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:51kDa.

## Clonality

Polyclonal